

General Features of the Mini Rechargeable Hearing Aid, HearAid.

Full Medix puts at your fingertips the product HearAid, is a mini hearing device type ITE, Inside The Ear. It is lightweight, comfortable, rechargeable and invisible to wear for a long time without pain.

Use the GM-916 hearing aid for older people who have hearing problems as long as they are not completely deaf. It can be used for near-deaf people in hearing aids with volume adjustment for nearsighted people. Deaf people do not need an ear test. Very easy to use.

The sound is constantly received and expanded, so there is the possibility of a howl after the ignition.

ENTRADA de carga: AC100-240V, 50-60Hz / 0.15^a

SALIDA del cargador: 5V, 1000mA

Chip: Analog chip-B

Max Output (DB): ≤110dB

Max Gain (DB): 30dB±5dB

Input Noise (DB): ≤15dB+3dB

Frequency Range (Hz): 300Hz-4500Hz

Total Harmonic Distort: ≤15%+3%

Battery Current: ≤4mA

Battery Type: Bulit in Recharge battery

Battery working standard time: About 18-20 Hours

Battery charge full: 1.5-2Hours

Sound adjustable: 5 levels

Material: ABS

Color: Baige

Size: G.W. : 2.6g, N.W. : 150g, Color box size : 12*9*4cm,

Carton size: 50.5*33*23.5cm

Packing: 1pcs

Certification: CE, FDA

Accessories: item*1pcs, earplug*5pcs, user manual*1, USB charger*1pcs, Brush*1

Applications of Mini Rechargeable Hearing Aid, HearAid.

HearAid Is suitable for the elderly and people with hearing loss.

Mini Rechargeable Hearing Aid Guarantee, HearAid.

HearAid has 1 year warranty subject to VentDepot clauses.



Características Técnicas Específicas del Mini Aparato Auditivo Recargable, HearAid.

Item #	Color	Frequency Range	Max Output	Voltage	Input Noise	Power		Material	Weight and Dimensions with Package in cm			
						Battery			kg	Wide	High	Deep
						Cable	Hearing Aid					
USFMA-001	Beige	300 to 4500Hz	110dB	1.2	≤15dB+3dB	Micro USB	7,50MAh	ABS Plastic	0.150	2.8	1.6	1.3

Specific Dimensions of the Mini Rechargeable Hearing Aid, HearAid in mm.

Item #	A	B	C
USFMA-001	28	16	13




Image gallery of the Mini Rechargeable Hearing Aid, HearAid.



ID	Author	Description Error	Description Update	Supervisor Name
V1	Angelica Romero	New Data Sheet	New Data Sheet	Jessica Lorenzo

