

## General Features of the Mini Rechargeable Hearing Aid, HearingMaxi

Full Medix puts at your fingertips the product HearingUltra, is a mini hearing device type ITE, Inside The Ear. It is lightweight, comfortable, rechargeable and invisible to wear for a long time without pain.

Anti-shock, fall and easy to use design.

Made from ABC of good quality, safe and non-toxic.

It can greatly reduce noise and provide clearer sound quality.

Suitable for seniors and people with hearing loss, it helps to hear the outside sound again, and regain confidence.

It contains a recharge case, to recharge the device anywhere and anytime.

With 300 mAh charging box, portable travel, charging anytime, anywhere. The headset can be charged 4-5 times

7,50MAh for single device. Type-C charge, can rechargeable, don't with power adapter.

Maximum OSPL90:  $\leq 118\text{dB} + 3\text{dB}$

High-frequency average SOPL90:  $103\text{dB} + 4\text{dB}$

Max Output (DB):  $\leq 110\text{dB}$

Max Gain (DB):  $30\text{dB} \pm 5\text{dB}$

Input Noise (DB):  $\leq 30\text{dB}$

Frequency Range (Hz):  $200\text{Hz} - 4000\text{Hz}$

Total Harmonic Distort:  $\leq 10\%$

Battery charge box: Charge product 4-5 times

Battery Type: Built in Recharge lithium battery

Working standard time: About 35Hours

Battery charge full: 2Hours

Sound adjustable: Smart touch

Material: ABS

Color: Black&White

Size: N.W.:3.3g, G.W. :140g, Color box size :  $12.5 \times 9.5 \times 4\text{cm}$ ,

Carton size:  $51.5 \times 48 \times 23.2\text{cm}$

Packing: 2pcs/box, 100pcs/carton, 14.2KG/CTN

Certification: CE, FDA

Accessories: item\*2pcs, charge case\*1, earplug\*10pcs, user manual\*1, USB cable\*1pcs, Brush\*1, etc.

## Applications of Mini Rechargeable Hearing Aid, HearingMaxi.

HearingMaxi is suitable for the elderly and people with hearing loss.

## Mini Rechargeable Hearing Aid Guarantee, HearingMaxi.

HearingMaxi has 1 year warranty subject to VentDepot clauses.



Características Técnicas Especificas del Mini Aparato Auditivo Recargable, HearingMaxi.

| Item #    | Color | Frequency Range | Max Output | Voltage | Input Noise | Power         |             | Material    | Weight and Dimensions with Package in cm |      |      |      |
|-----------|-------|-----------------|------------|---------|-------------|---------------|-------------|-------------|--|------|------|------|
|           |       |                 |            |         |             | Battery       |             |             | kg                                       | Wide | High | Deep |
|           |       |                 |            |         |             | Cable         | Hearing Aid |             |  |      |      |      |
| MXFMH-001 | White | 200 to 4000Hz   | 110dB      | 1.2     | ≤30dB       | Type-C charge | 7,50MAh     | ABS Plastic | 0.200                                    | 51.5 | 48   | 23.2 |
| MXFMH-002 | Black | 200 to 4000Hz   | 110dB      | 1.2     | ≤30dB       | Type-C charge | 7,50MAh     | ABS Plastic | 0.200                                    | 51.5 | 48   | 23.2 |

Specific Dimensions of the Mini Rechargeable Hearing Aid, HearingMaxi in mm.

| Item #    | A    | B  | C  | D   | E  | F  |
|-----------|------|----|----|-----|----|----|
| MXFMH-001 | 11.5 | 20 | 10 | 125 | 95 | 40 |
| MXFMH-002 | 11.5 | 20 | 10 | 125 | 95 | 40 |





Image gallery of the Mini Rechargeable Hearing Aid, HearingMaxi.



USB rechargeable, used in various scenarios



Power Bank



Charger



Computer



## AUTOMATIC NOISE REDUCTION


Clear sound, effectively reduce noise. Cut off howling, smarter sound processing



Touch control operation, simple and convenient



Image gallery of the Mini Rechargeable Hearing Aid, HearingMaxi.



| ID | Author          | Description Error | Description Update | Supervisor Name |
|----|-----------------|-------------------|--------------------|-----------------|
| V1 | Angelica Romero | New Data Sheet    | New Data Sheet     | Jessica Lorenzo |