

General Features of the Mini Rechargeable Hearing Aid, LoudAudi.

Full Medix puts at your fingertips the product LoudAudi, is a mini hearing device type, Inside the Ear. It is lightweight, comfortable, rechargeable, and invisible to wear for a long time without pain.

Use the GM-916 hearing aid for older people who have hearing problems as long as they are not completely deaf. It can be used for near-deaf people in hearing aids with volume adjustment for nearsighted people. Deaf people do not need an ear test. Very easy to use.

The sound is constantly received and expanded, so there is the possibility of a howl after the ignition.

Maximum OSPL90: ≤110dB

High-frequency average SOPL90: 106dB

Max Output (DB): ≤110dB

Max Gain (DB): 30dB±5dB

Input Noise (DB): ≤25dB

Frequency Range (Hz): 200Hz-5000Hz

Total Harmonic Distort: ≤10%

Battery charge box: Charge product 4-5 times

Battery Type: Built in Recharge lithium battery

Working standard time: About 30 Hours

Battery charge full: 2Hours

Material: ABS

Color: Red&Blue

Size: N.W.:4.5g, G.W.:160g, Color box size : 12.5*9.5*4cm,

Carton size: 51.5*48*23.2cm

Packing: 2pcs

Accessories: item*2pcs, charge case*1, user manual*1, USB cable*1pcs, Brush*1, etc.



Applications of Mini Rechargeable Hearing Aid, LoudAudi.

HearAid Is suitable for the elderly and people with hearing loss.

Mini Rechargeable Hearing Aid Guarantee, LoudAudi.

HearAid has 1 year warranty subject to VentDepot clauses.

Características Técnicas Específicas del Mini Aparato Auditivo Recargable, LoudAudi.

Item #	Color	Frequency Range	Max Output	Voltage	Input Noise	Power		Material	Weight and Dimensions with Package in cm			
						Battery	Cable		kg	Wide	High	Deep
						USLDA-001	Red&Blue		200Hz-5000Hz	110dB	1.3	≤25dB



Specific Dimensions of the Mini Rechargeable Hearing Aid, LoudAudi in mm.			
Item #	A	B	C
USLDA-001	120	90	40




Image gallery of the Mini Rechargeable Hearing Aid, LoudAudi.



MAGNETIC CHARGING LONG BATTERY LIFE

Type-c charging is faster and can be inserted regardless of the positive or negative







-  **2-3H**
Charing time
-  **30H**
Working time

When charging, the display will show the charging status and battery status




PRODUCT ADVANTAGES

This is a high-end fashion hearing aid

-  Digital Display
-  Comfortable
-  Lightweight
-  Intelligent
-  Noise reduction
-  Fashion




Automatic noise reduction sound quality is more delicate and clear



Push button design one key operation, the elderly simple to understand



On/off key
Volume knob



ID	Author	Description Error	Description Update	Supervisor Name
V1	Angelica Romero	New Data Sheet	New Data Sheet	Jessica Lorenzo

